

## SUPERPOLO 220



### Characteristics

<b>Temperature Range</b>	-10 /-45 °C <i>Also available PLUS version -10 / -60°C</i>
<b>Capacity</b>	230 lt (8,13 cu ft)
<b>External Dimensions (W x D x H)</b>	cm 105,5 X 73 X 86,5 <i>Internal Dimensions (W x D x H) cm 85x45x65 (+ 20x45x40); Weight kg 64</i>
<b>Structure and Insulation</b>	embossed aluminium OR white coated steel inside and sheet steel epoxy coated white colour outside, thickness 100 mm
<b>Feet</b>	n.4 feet
<b>Door</b>	nr 1 lid
<b>Key Lock</b>	
<b>Controller</b>	Digital electronic thermostat
<b>Alarms</b>	Acoustic and Visual for temperature deviations
<b>Refrigeration Type</b>	Static
<b>Defrost</b>	Manual

### Technical data

<b>Details</b>	<ul style="list-style-type: none"> <li>• Refrigerant type: R290 CFC free</li> <li>• Mains Voltage 230V/1/50Hz</li> <li>• Power consumption 24H: 5,1 kWh/24H</li> <li>• Instant power consumption: 0,22 kW (2,2 A)</li> <li>• Noise level &lt; 50 dB</li> <li>• Ambient Temperature: +10/+32 °C; humidity 60%</li> <li>• Connector: Schuko</li> </ul>
<b>In accordance with</b>	2014/30/CE, 2014/35/CE and following modifications
<b>Package</b>	Strong Cardboard on wooden pallet <i>Dim. cm.110x80x110; Kg.74</i> <ul style="list-style-type: none"> <li>• Wooden crate on wooden pallet available: dim. cm.125x90x112; Kg.106</li> </ul>

### Accessories

Castors (2 with with brakes)

DMLP UNIT

*External module with DMLP DIGITAL MONITOR: Monitor/alarm/recording system with rechargeable battery*  
*- Temperatures and alarms recorder on removable Secure Digital (SD Back-Up memory): it can store 5 years data, recording temperature and alarms every minute. The SD card can be connected to any PC: a software (Capture) is included, so that you can download, save and print all data from your PC*  
*- Battery backup: ensures alarms and recording of temperatures even in case of power failure. The battery system lasts about 12 hours.*  
*- Dry contact for remote connection + RS485 port*

DMLP PRINTER

Internal steel wire basket

*dimensions W x D x H cm 40x25x15*

Voltage stabilizer

60 Hz Frequency

115 V

**FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy

Tel. +39 0522 976232- Fax +39 0522 976028

info@fiocchetti.it - www.fiocchetti.it

P.I.00934960352

Registro A.E.E. IT8020000002833