

B65W



Characteristics

External Dimensions (W x D x H)	cm 53,5 x 60 x 91 / 103,3
Structure and Insulation	stainless steel inside and outside
Height adjustable Feet	
Kind of ice	cube
Motor cooling	water
Production Kg/24h (up to)	75
Bin capacity (kg)	40
Weight	kg 45

Technical data

Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption kWh/24h: 13 kW/h
- Water consumption l/h: 41
- Fuse (A): 10
- Cooling Capacity Btu/h/W: 1570
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 58 dB
- Ambient Temperature: +18/+43 °C; humidity 60%

In accordance with	2014/30/CE, 2014/35/CE and following modifications
Package	Strong Cardboard on wooden pallet <i>Dim. cm.76 x 63 x 113; Kg.52</i> <ul style="list-style-type: none">• Wooden crate on wooden pallet available: dim. cm.87 x 69 x 112; Kg. 68
note	Available also air cooled

Accessories

60 Hz Frequency

Wooden crate Package

FIOCCHETTI SCIENTIFIC S.R.L.

Via Panagulis, 48 - 42045 Luzzara (RE) Italy
Tel. +39 0522 976232- Fax +39 0522 976028
info@fiocchetti.it - www.fiocchetti.it
P.I.00934960352
Registro A.E.E. IT8020000002833