

## B90A



### Characteristics

<b>External Dimensions (W x D x H)</b>	cm 70 x 60 x 97 / 109
<b>Structure and Insulation</b>	stainless steel inside and outside
<b>Height adjustable Feet</b>	
<b>Kind of ice</b>	cube
<b>Motor cooling</b>	air
<b>Production Kg/24h (up to)</b>	93
<b>Bin capacity (kg)</b>	50
<b>Weight</b>	kg 58

### Technical data

#### Details

- Refrigerant type: R134A CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption kWh/24H: 18
- Water consumption l/h: 9
- Fuse (A): 16
- Cooling Capacity Btu/h/W: 7745
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 58 dB
- Ambient Temperature: +10/+32 °C; humidity 60%

<b>In accordance with</b>	2014/30/CE, 2014/35/CE and following modifications
<b>Package</b>	Strong Cardboard on wooden pallet  <i>Dim. cm.76 x 63 x 106; Kg.68</i>  • <i>Wooden crate on wooden pallet available: dim.dim. cm.86 x 68 x 111; Kg. 77</i>
<b>note</b>	Available also water cooled

### Accessories

60 Hz Frequency

Wooden crate Package

#### **FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy  
Tel. +39 0522 976232- Fax +39 0522 976028  
info@fiocchetti.it - www.fiocchetti.it  
P.I.00934960352  
Registro A.E.E. IT8020000002833