

## BF80W



### Characteristics

**External Dimensions (W x D x H)** cm 53 x 63 x 92,1

**Structure and Insulation** stainless steel inside and outside

**Height adjustable Feet**

**Kind of ice** flakes

**Motor cooling** water

**Production Kg/24h (up to)** 73

**Bin capacity (kg)** 25

**Weight** kg 51

### Technical data

#### Details

- Refrigerant type: R134a CFC free
- Mains Voltage: 230V/1/50Hz
- Power Consumption: 7,5 kW/h
- Water consumption l/h: 16
- Electrical consumption: 330 W
- Fuse (A): 10
- Cooling Capacity Btuh/W: 670
- Water temperature min/max: 5°C-35°C
- Pressure min/max: 1 bar/5 bar
- Noise level < 60dB
- Ambient Temperature: +18/+43 °C; humidity 60%

**In accordance with** 2014/30/CE, 2014/35/CE and following modifications

**Package** Strong Cardboard on wooden pallet

*Dim. cm.60 x 69 x 98; Kg.60*

- *Wooden crate on wooden pallet available: dim. cm.66 x 79 x 102; Kg. 80*

**note** Available also air cooled

### Accessories

60 Hz Frequency

Wooden crate Package

#### **FIOCCHETTI SCIENTIFIC S.R.L.**

Via Panagulis, 48 - 42045 Luzzara (RE) Italy  
Tel. +39 0522 976232- Fax +39 0522 976028  
info@fiocchetti.it - www.fiocchetti.it  
P.I.00934960352  
Registro A.E.E. IT8020000002833